SOGO

Human Technology

manual de instrucciones / instructions for
use mode d'emploi / manual de instruções
gebrauchsanweisung / manuale d'istruzioni

Almohadilla Eléctrica
Electric heating Pad
Oreiller électrique
Almofada Elétrica
Heizkissen
Pad elettrico

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Instructions for use

PLEASE READ ALL THE INSTRUCTIONS PRIOR TO USING YOUR HEATING BLANKET.

SAFETY GUIDELINES

This manual can also be downloaded from our website: www.sogo.es

1. Do not fold the heating pad. Do not use the pad in folded form.
2. Do not use pins on the heating pad.
3. Do not wrinkle the heating pad.
4. Do not puncture, cut or otherwise damage any part of the heated electric pad.
5. Keep away from water, moisture and high humidity.
6. Do not use when the pad is wet or moist. Use only when it is completely dry.
7. Do not use with a helpless or insensitive person or infant.
8. Do not attempt to modify, repair, disassemble or reassemble any component(s) of the heated electric pad. Repairs should be carried out by the manufacturer or other qualified persons.
9. This appliance is not intended for use by persons (including children) whose physical, sensory or mental capabilities are impaired, or by persons lacking experience or knowledge, unless they have had the benefit of prior supervision and instruction in using the device from a person responsible for their safety. Children should be supervised to ensure that they do not play with the appliance.
10. This product is not intended for use in

hospitals. It is for domestic use only.

11. Warning! Never go to sleep with the heating pad in use. Extended use can result in skin burns.

12. Do not cut or pierce the heating pad.

13. Do not pull, twist or otherwise apply force to the cables.

14. The device should not be connected to the mains via a timer.

15. The switch should not be placed on or under the appliance, nor covered when the appliance is in operation.

16. When in use do not put heavy or large objects on the pad

17. When in use do not place any source of heat over the pad, e.g. hot water bottles, heating pads or similar

18. When not in use, unplug the pad from the socket. The appliance should be stored in a dry place. Take care not to put additional weight on it.

19. Do not allow children being supervised to play with the appliance.

20. If the supply cord is damaged, it must be replaced by the manufacturer, its service agent or similarly qualified persons in order to avoid a hazard.

21. The appliance is not to be used if there are signs of damage. Observe for signs of damage such as rips, stitching coming undone, cords becoming loose or similar.

22. The appliance is not to be used for warming animals.

23. The appliance has to be repaired or replaced if the cover is worn.

PRODUCT FEATURES AND TECHNICAL SPECIFICATIONS

For relieving pain in the neck, shoulders and the back

Particularly indicated for osteoarthritis and rheumatism treatments

Material: Super Polar fleece cover with inner polyester + PVC Heating Pad

With electronic temperature control with two settings of High/Low

LED indicator on switch control

Detachable and machine washable cover

With overheating protection system

Size: 30 x 40cm

Voltage: 230V~50/60HZ

Power: 45W

Net Weight 400 gms

PARTS OF THE CONTROLLER

- 1) 2 Levels of heat to select from
- 2) LED indicator

OPERATION

- Plug into the mains socket and switch on at the plug.
- Adjust the temperature setting to the desired level:
0 = Off (disconnection)
1 = Low heat
2 = High heat
- The heated Pad / Heating Jacket takes about 10 minutes to reach the desired temperature and then automatically regulates its operating temperature.
- Safety cut off temperature of 60 °C.

INSPECTION AND MAINTENANCE

- The heated electric pad should be checked regularly to ensure the control box and the pad is intact. Use of the pad should be stopped immediately if any wiring is exposed or the pad / control box is compromised.
- **The heated electric pad should be kept in a cool dry place.**
- **DO NOT EXPOSE THE CONTROLLER, CABLES OR HEATING PAD TO WATER, MOISTURE OR HIGH HUMIDITY**
- Sponge clean only with luke-warm water and neutral detergent. Rinse with clean water by sponging and allow to air dry before re-use.
- **WARNING:** do not allow the electric control to get wet while cleaning.

Small stains can be cleaned with a cloth or a damp sponge, and possibly with a small amount of washing liquid for delicate fabrics. **Caution:** do not dry clean the jacket, wring it out, dry it in a tumble dryer or a mangle and do not iron.

The cover of the heating pad can be washed in a washing machine. Use a cycle for delicate fabrics at 30°C. Use a washing product for delicate fabrics and use the amount recommended by the manufacturer. **Do not wash the jacket too often as this could cause damage. It should not be washed in a washing machine more than 5 times during its lifetime.**

To dry, do not hang up using clothing pegs or any other items.

Do not expose to sources of heat such as the sun, radiators etc.

Note.: Under no circumstances you should not switch on the jacket to dry it.

STORAGE

If you are not going to use the heating pad for a long period, we recommend storing it in its original packaging in a dry place and that you do not place anything on top of it.

Before storing the heating pad make sure to follow the following steps:

- Make sure the control is in «position 0».
- Disconnect the power supply.
- Leave the heated pad to cool down completely.
- Store heated pad preferably without folding with the cable and plug folded carefully without damaging it and keep it in a protected place.

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DECLARATION OF CONFORMITY

This device complies with EU Low Voltage Directive 2014/35/EC.
Electromagnetic Compatibility Directive 2014/30/EU.
Directive 2015/863/EU on the restriction of the use of certain hazardous substances in electrical.
Directive 2009/125/EC on the eco-design requirements applicable to energy-related products.

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This symbol on the product or on the packaging indicates that this product can't be disposed as normal rubbish or household waste. All the electrical, electronic equipment's and battery-operated units must recycle in proper manner and according to the local municipal laws. You can recycle them by taking them to government authorized disposal centres or specialized bins which you can find in any nearby big super markets, electronics or electro domestics products stores or malls who have these types of facilities available.

Designed by: SOGO based on European quality standards
Imported by: Sanyan Appliances S.L, NIF: B98753056, C/ Barcas 2, 2, 46002 Valencia, Spain
Product manufactured in CHINA. After-sales service: www.sogosat.com sogosat@sogosat.com / 0034 902 222 161

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