SOGO®

Human Technology

ES Manual de instrucciones
EN Instructions for use
FR Mode d'emploi
PT Manual de instruções
DE Gebrauchsanweisung
IT Manuale d'istruzioni
CA Manual d'istruccions
CZ Návod k použití

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1. IMPORTANT NOTICE

- Always read the instruction book carefully before using.
- This manual can be downloaded from our web page www.sogo.es
- Keep these instructions for future reference.

2. SAFETY INSTRUCTIONS FOR THE USER

GENERAL PRECAUTIONS

GENERAL PRECAUTIONS

I. Restrictions

- Never immerse the device or the charging cable in water or any other liquid. The unit is not waterproof. In the event of the appliance falling in water take it to an authorized service agent for repair before reusing.
- Do not place or operate the device on or near hot surfaces (e.g., stove plates) or open flames.
- In case of hardware problems, do not attempt to repair the product yourself. Repairs should only be carried out by qualified technicians.
- Never use accessories that are not recommended by the manufacturer.
- The use of accessories not recommended or sold by the appliance manufacturer may result in fire, electric shock or irreparable damage to the product.
- Do not use the device for any other purpose than described in this manual.
- Do not use the device while in charge, otherwise it may also damage the device.
- Do not drop, disassemble, deform, modify or insert any objects into this product as operating damage may result.
- Do not attempt to repair, modify or disassemble the product yourself. This product does not contain any user-serviceable components and may void the warranty.
- Do not clean the product with alcohol, ammonia-based cleaners and abrasive cleaners.
- Do not place this product in to contact with any sharp objects as this may cause scratches and operating damage.

- Do not attempt to replace the product battery, it is built-in and is not user changeable.
- Cleaning and user maintenance shall not be made by children.
- This Product is not a toy, keep it out of reach of children.
- Keep this product out of reach of children and pets to avoid chewing and swallowing small parts.
- This appliance is not intended for use by persons (including children) with reduced physical sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety. Children should be supervised to ensure that they do not play with the appliance.

II. Precautions while using this device

- Keep the device away from direct sunlight, fire, or other high voltage equipment.
- **WARNING:** the device might get warm during charging.
- Do not replace the attachment head while it is running.
- Avoid exposure to extreme temperature, direct sunlight, moisture, sand, dust, or mechanical shock; never leave the products directly into the sun.

Medical Warning

This massage gun is intended for muscle relaxation only. It is not a medical device and is not intended to diagnose, treat, cure or prevent any disease.

Consult your physician before use if you:

- Are pregnant.
- Have diabetes or neuropathy.
- Wear a pacemaker or other implanted medical device.
- Have recently undergone surgery.
- Have any medical condition that may be affected by vibration massage.

Do not use the massage gun:

- On the head, chest or genital area.
- On wounds, swollen or inflamed areas.

- Directly over the spine or joints.

Stop using immediately if pain or discomfort occurs.

Battery & Charging Safety

- Use only a 5V/2A USB charger or a compatible fast charger.
- Fully charge the battery before first use.
- The device may become warm during charging.
- Recharge the battery promptly when the low battery indicator appears.
- Do not dispose of batteries in fire.

3. TECHNICAL SPECIFICATIONS

Power: 42 W

Speed settings: 5

Maximum speed: 3000 rpm

Maximum clamping force: 9 kg

Deep-tissue massage amplitude: 6 mm

Automatic shut-off: after 10 minutes of use

Charging time: 3.5 hours

Operating time: 4.5 hours max. (minimum speed)

LED display

Touch-sensitive control button

3.7 V, 1800 mAh lithium battery

Type-C charging port

Charging indicator: flashing white light

Battery indicator: 5% < battery < 20%, flashing red light

Low battery indicator: < 5%, steady red light

Noise level: 60 dB

Charging Adaptor: 5 V/2 A (not included)

Accessories: Carrying case, charging cable, 3 massage heads (round silicone head, U-shaped head and flat plastic head)

This product complies with EU regulation 2023/826 under directive 2009/125/EC on the Eco design requirements.

Power consumption in off mode: 0.0W

Power consumption in standby mode: N/A

Time after which the product enters into off mode: N/A

4. PARTS DESCRIPTION

IMAGE 1

MASSAGER

1. Massage head socket
2. Handle
3. speed indicator lights
4. Power / Speed button
5. Battery indicator light
6. Charging port

Accessories

Attachment heads

IMAGE 2

1. Ball head
2. Fork head
3. Flat head
4. USB charging cable
5. Carrying case

5. INSTRUCTIONS FOR USE

BEFORE USE

1. CHARGING INSTRUCTIONS

- Connect the USB Type-C charging cable to the massage gun and then connect it to a 5V/2A USB charger (not included).
- The battery requires approximately 3.5 hours to fully charge.
- The massage gun cannot be operated while charging.
- Recharge the battery whenever the low battery indicator is displayed.

2. Installing the Massage Head

- Choose the massage head suitable for the desired muscle group.
- Firmly insert the massage head into the socket until fully seated.
- Always switch off the massage gun before replacing any attachment.
- Never change the massage head while the device is operating.

3. Power On/Off

Press and hold the Power Button for approximately 1.5 seconds to switch the massage gun ON or OFF.

4. Speed selection

While the massage gun is operating, briefly press the Power Button to cycle through the available massage levels.

The massage gun automatically remembers the last speed level used and starts at the same level the next time it is switched on.

Levels	Indicator light status	
Auto	4 white lights	Smart pressure control
Level 1	1 white light illuminates	Muscle awakening
Level 2	2 white lights illuminates	Fascia relaxation
Level 3	3 white lights illuminates	Activation for exercise
Level 4	4 white lights illuminates	Deep tissue massage

5. Memory mode

The device intelligently remembers the last speed/level you used and automatically enters in that speed while turning on next time.

6. PROTECTION

Over-pressure Protection

If excessive pressure is applied, the massage gun will temporarily stop operating to protect the motor. Reduce the pressure and press the Power Button to continue.

Auto Shut-off

For safety, the massage gun automatically switches off after 10 minutes of continuous use. Simply switch it on again if further massage is required.

7. INDICATOR LIGHTS

Status	Battery indicator status	Meaning
Battery	Flashing red light	5% < Battery Level < 20%
	Steady red light	Battery Level ≤ 5% (Low)
Charging	Flashing white light	Charging
	Steady white light	Fully charged

8. MASSAGE AREAS

Area	Time
	Back: 10 ~ 20s x 3 times
	Arm: 10 ~ 20s x 3 times
	Buttocks: 20 ~ 60s x 4 times
	Legs: 10 ~ 45s x 2 times
	Soles of the feet: 20 ~ 30s x 3 times

NOTE:**Tips:**

- Do not massage directly over the spine.
- Avoid bones, joints and sensitive areas.
- Do not massage wounds or swollen muscles.
- Do not use behind the knees, inside the armpits or near the collarbone.
- Massage slowly following the direction of the muscle.
- If you experience pain or discomfort, stop using immediately.

6. CLEANING AND MAINTENANCE

- Switch off the massage gun before cleaning.
- Remove the massage head before cleaning.
- Wipe the appliance with a soft, dry or slightly damp cloth.
- Do not immerse the massage gun in water.
- Do not use alcohol, solvents or abrasive cleaners.
- Store the product in a cool, dry place away from direct sunlight.
- Keep the carrying case clean and dry.

7. DISPOSAL OF ELECTRONIC GOODS

DISPOSAL OF BATTERY-OPERATED UNITS

Do not dispose the battery-operated units as normal rubbish, always these should be recycled in proper manner and according to the local municipal laws. You can recycle them by taking them to government authorized disposal centres or specialized bins which you can find in any nearby big super markets, electronics or electro domestics products stores or malls who have these types of facilities available.

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This symbol on the product or on the packaging indicates that this product can't be disposed as normal rubbish or household waste. All the electrical, electronic equipment's and battery-operated units must recycle in proper manner and according to the local municipal laws. You can recycle them by taking them to government authorized disposal centres or specialized bins which you can find in any nearby big super markets, electronics or electro domestics products stores or malls who have these types of facilities available.

Designed by: SOGO based on European quality standards

Imported by: Sanytan Appliances S.L, NIF: B98753056, C/ Barcas 2, 2, 46002 Valencia, Spain

Product manufactured in CHINA. After-sales service: www.sogosat.com sogosat@sogosat.com / 0034 902 222 161

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