SOGO®

Human Technology

ES Manual de instrucciones
EN Instructions for use
FR Mode d'emploi
PT Manual de instruções
DE Gebrauchsanweisung
IT Manuale d'istruzioni
CA Manual d'instruccions
CZ Návod k použití

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1. IMPORTANT

- Always read the instruction book carefully before using.
- This manual can be downloaded from our web page www.sogo.es
- Keep these instructions for future reference.

2. SAFETY INSTRUCTIONS FOR USER

- Place the Bathroom scale on flat, hard and stable floor surface.
- The max. weighing limit of the scale is 180 kg, do not exceed this limit.
- To ensure the life of the bathroom scale, do not keep it under wet, hot & cold condition.
- The product belongs to high precision electronic weighing apparatus, please do not press or drop the scale otherwise the scale is easy to damage.
- If the scale surface is dirty, please clean the scale with soft cloth, don't use any chemical detergent or allow water to get inside the scale.
- This product is for personal body analysis only. Do not use for medical diagnosis.
- If you think the results are incorrect, then check the batteries, if they are connected well with both of poles and in case not, place them correctly. Still if you don't get the correct results, then change the batteries and try again.
- Not to use for the persons with pacemaker.
- Do not use on slippery floors
- Do not use when your body is wet.
- Do not use it with rechargeable 1.2V, AAA Batteries.
- People with implanted medical devices (such as pacemakers or defibrillators) must not use the body composition measurement function.
- Not suitable for persons with limb amputations or severe foot contact limitations.
- Do not jump on the scale.
- It is recommended to measure at the same time each day under the same conditions.

3. PRODUCT FEATURES

High resistance Tempered glass base

LED Display

Multi frequency BIA technology

Bluetooth connected (FEELFIT APP)

Support iPhone & Android

Display size: 11.6 x 7.2cm

Electrodes: 8

Capacity: 0.2 kg to 180 kg

Graduation: 100g

Selectable Unit: KG/LB/ST

Key body parameters: 8

App parameters 50

Auto Power On/Off

Overload & Low Battery Indication

Power: 3 x AAA (Not Included)

Dimensions: 30 x 32.4 x 3.1 cm

Measures Body fat, Water, Muscle Mass, Bone Mass, Visceral fat, BMI etc

How to Ensure Accurate Readings

- Always weigh barefoot. Do not step on the scale with wet feet, wearing socks or shoes to prevent slipping, damage and inaccurate readings.
- Recommended to be placed on hard and flat surface only.
- Always weigh yourself on the same scale, same place and same surface if possible.
- Recommended measurement time should be after waking up or before the breakfast and next weight should be two hours after your meal but before your next meal, or two hours after dinner but before bed.
- Wrong measurements can occur if feet are dirty, dry or extremely callous.
- Body fat percentage estimates will vary with the amount of water in the body and can be affected by dehydration or over-hydration due to some factors as alcohol consumption, menstruation, illness, intense exercise.
- Stand upright and not to the edge of the scale to ensure good contact with the film or electrodes (the film and electrodes are both conductive media depending on the scale models).
- Horizontal line should be running below arch of each foot.
- Complete features are available only with the apps.

4. PARTS DESCRIPTION

A. Top view

1. Thumb electrodes
2. Digital display
3. Foot electrodes area

B. Bottom view

1. Hand electrodes
2. Handle cable
3. Handle cable port
4. Battery compartment

5. INSTRUCTIONS FOR USE

INITIAL SETUP

Battery installation

Open the battery compartment located on the underside of the scale and insert 3 x AAA batteries, ensuring the correct polarity (+/-). Close the battery cover securely. When the display shows "Lo", replace all batteries with new ones.

Step 1: Download app and create an account

Scan the QR code provided or download the FEELFIT app from the App Store or Google Play. Install the app and register a user account by following the on-screen instructions.

Step 2: Pair scale with smartphone.

Enable Bluetooth® on your smartphone and open the FEELFIT app. Select Add Device "SS-2950" and follow the instructions to pair the scale with your phone. Once paired, the scale will automatically reconnect during future measurements when Bluetooth is enabled.

How to use

Stepping onto the Body Fat Scale

Place the scale on a hard, flat surface. Step onto the scale barefoot while holding the handle with both hands. Remain still until the measurement is complete. Your body composition data will be displayed on the scale and synchronised with the FEELFIT app.

Note: The scale can be used without the app to measure body weight only. To view body composition data and measurement history, the scale must be connected to the FEELFIT app.

Body fat scale measurement demonstration, **BODY POSITION**

For accurate body composition analysis, stand upright in a relaxed position while holding the handle firmly with both hands. Keep your arms slightly away from your body and avoid touching your waist or thighs. Remain still and maintain the same posture until the measurement is completed.

Important:

- Keep your arms slightly extended away from your body.
- Do not bend your elbows or wrists.
- Do not move or talk during the measurement.
- Wear light, close-fitting clothing whenever possible.

Body fat scale measurement demonstration, **FEET POSITION**

Correct foot placement is essential to ensure proper contact with the electrodes and accurate measurement results. Always stand barefoot with clean and dry feet, placing both feet fully on the electrode areas.

Important:

- Keep both feet parallel.
- Ensure the entire sole is in contact with the electrodes.
- Keep your heels separated.
- Avoid touching your inner thighs together.
- Do not stand on the edges of the scale.

Body fat scale measurement demonstration, **HANDS POSITION**

Hold the handle comfortably with both hands, ensuring full contact between your fingers, thumbs and the electrodes. Proper hand placement allows the electrical signal to pass correctly through the upper body for a complete body composition analysis.

Important:

- Place your thumbs on the round electrodes.
- Place your fingers on the long electrodes.
- Keep the handle cable facing downward.
- Do not clench your hands.
- Do not release or reposition the handle until the measurement is complete.

Skin Contact

For reliable body composition measurements, ensure that your hands and feet are clean, dry and free from lotions or creams. Poor skin contact may result in unstable or inaccurate readings.

For best results:

- Measure barefoot.
- Clean and dry your hands and feet before use.
- Avoid measuring immediately after intense exercise, bathing or using hand cream.
- Take measurements at approximately the same time each day under similar conditions.

Body metrics display

The scale displays weight, body fat, body water, muscle mass, bone mass, visceral fat, BMI, and rotation (switching between weight and BMI data, each displayed for 3 seconds). For additional metrics, please refer to the app after measurement.

Icon	Details
	Body fat percentage icon: illuminates when displaying body fat percentage measurement results.
	Body water icon: illuminates when displaying body water measurement results.
	Muscle icon: illuminates when displaying muscle mass measurement results.
	Bone mass icon: illuminates when displaying bone mass measurement results.
	Visceral fat level icon: illuminates when displaying visceral fat level measurement results.
	BMI icon: illuminates when displaying BMI measurement results.
	Percentage unit

Display Messages

Display	Meaning
Lo	Low Battery
Err	Overload
Err1	Unstable Measurement
----	Try Again
No Display	Replace Batteries
Bluetooth Icon	When it flashes, check the device, icon steady, Bluetooth connected

6. CLEANING AND MAINTENANCE

- Wipe with a soft damp cloth.
- Do not immerse in water.
- Keep dry when not in use.
- Store in a cool, dry place.
- Remove batteries if not used for a long period.

7. BATTERY DISPOSAL**DISPOSAL OF BATTERY-OPERATED UNITS**

Do not dispose the battery-operated units as normal rubbish, always these should be recycled in proper manner and according to the local municipal laws. You can recycle them by taking them to government authorized disposal centres or specialized bins which you can find in any nearby big super markets, electronics or electro domestics products stores or malls who have these types of facilities available.

DISCLAIMER:

This product is compatible with the FEELFIT app, which is developed and maintained independently by a third-party provider. SOGO is not affiliated with the app developer and accepts no responsibility for the app's availability, functionality, updates, data storage, privacy practices or future compatibility. The app's terms of use and privacy policy are governed solely by the app provider.

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Alimentation : 3 piles AAA (non fournies)

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This symbol on the product or on the packaging indicates that this product can't be disposed as normal rubbish or household waste. All the electrical, electronic equipment's and battery-operated units must recycle in proper manner and according to the local municipal laws. You can recycle them by taking them to government authorized disposal centres or specialized bins which you can find in any nearby big super markets, electronics or electro domestic products stores or malls who have these types of facilities available.

Designed by: SOGO based on European quality standards

Imported by: Sanyan Appliances S.L, NIF: B98753056, C/ Barcas 2, 2, 46002 Valencia, Spain

Product manufactured in CHINA. After-sales service: www.sogosat.com sogosat@sogosat.com / 0034 902 222 161

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